



Ask About Web Special !
Call 800 788-2011

H/S Catheter Set

For enhanced visualization of uterine pathology

Features & Benefits

- Available with integrated stylet for easy insertion
- Unique placement sheath eliminates the need for a tenaculum
- Comes in both 5 and 7 French
- Ideal for office use
- Latex-free construction
- The Original Ackrad Catheter
- Over 30 Years as the Industry Leader
- Made in America 



Sonohysterography has been proven to significantly improve the depiction of uterine pathology

CooperSurgical

*Keeping You at the Forefront
of Women's Health Care™*

Ordering information

H/S Catheter

The Uterine Infusion Catheter family comes in a standard configuration or with an integrated stylet. Both versions are minimally invasive, easy-to-use devices that ensure a diagnostically enhanced*, yet gentle, patient exam. Designed for use in saline-infusion sonohysterography and hysterosalpingography.

H/S Catheter with Integrated Stylet

Based on the original Ackrad round, latex-free, balloon catheter design, this is used when additional catheter rigidity is needed.

- #61-3005 H/S Catheter with integrated Stylet - 5 French, 10 per box
- #61-3007 H/S Catheter with integrated Stylet - 7 French, 10 per box

Standard H/S Catheters

- #61-5005 H/S Catheter, Standard - 5 French, 10 per box
- #61-5007 H/S Catheter, Standard - 7 French, 10 per box

H/S Stylet™


- #61-0020 H/S Stylet, 10 per box

SIS is useful for finding focal lesions in AUB

Pre-Saline Infusion



Post-Saline Infusion



	Sensitivity	Specificity
Transvaginal Ultrasound	75%	76%
Saline Infused Sonohysterography	93%	94%

*Kamel, HS, Darwish, AM, Mohamed, SA. Comparison of transvaginal ultrasonography and vaginal sonohysterography in the detection of endometrial polyps. Acta Obstet Gynecol Scand 2000; 79:60



H/S Procedure Tray - Product No. 61-5205

Each Procedure Tray contains:

- Fenestrated drape
- Disposable speculum
- Disposable cervical dilator
- Three swab sticks
- One wide swab
- Sterile wrap
- Two 3" x 3" gauze pads
- Inner tray for Povidone-iodine solution
- Povidone-iodine ointment
- Two packs of Povidone-iodine solution
- 36" extension tube
- 20 cc syringe with 18-gauge needle
- Lubricating jelly packet
- Sanitary napkin

Complete H/S Procedure Trays

- H/S Procedure Tray - Round 5 French - Product number: 61-5205
- H/S Procedure Tray - Round 7 French - Product number: 61-5207
- H/S Procedure Tray - Round 5 French with Stylet - Product number: 61-3205
- H/S Procedure Tray - Round 7 French with Stylet - Product number: 61-3207

@operSurgical



Heartland Medical Specialties, Inc.
9854 Industrial Blvd
Lenexa, Kansas 66215
Fax: 913 541-1897
www.heartlandmedspec.com

To find out more about the H/S Catheter, or to place an order please call 1 800 788-2011

Form # 80972 03/08

© Copyright 2008 CooperSurgical, Inc