

Simple and effective way
to allow visibility of the
uterine cavity



FEATURES

- Placement sheath decreases the need for a tenaculum
- The distal end of the catheter is balloonless allowing for a more complete view of the uterus during procedure
- An ink band provides a reference mark to indicate that proper catheter positioning has been received
- The round, closed catheter tip eases the introduction of the catheter into the uterus
- Sliding sheath assist in catheter placement

 Not made with natural rubber latex. PVC Free & non-DEHP

TECHNICAL DATA

French Size	5Fr	Maximum Balloon Volume	5Fr, 1.5ml/7Fr, 3ml
Standard Length	28cm	Box Qty.	10 ea.
Sterilization Method	ETO		

ORDERING INFORMATION

TMI1170 Sono-Inject Catheter, 5Fr

Download IFU





Simple and effective way
to allow visibility of the
uterine cavity

FEATURES

- Specially designed tip helps prevent leakage
- Rigid catheter sheath allows for simple placement and single operator use
- Long tip at distal end improves seating in cervical canal and has a wide tip opening for increased flow
- Balloonless catheter allows for a more complete view of the uterine cavity during procedure

 Not made with natural rubber latex. PVC Free & non-DEHP

TECHNICAL DATA

French Size	9 Fr	Sterilization Method	ETO
Length	30 cm	Box Qty.	10 ea.

ORDERING INFORMATION

TMI1162-HF eZInject SIS Catheter with High-Flow




Catheter Tip Length



Catheter Tip Opening



Heartland Medical Specialties, Inc.
8347 Melrose Drive
Lenexa, KS 66214
Telephone: 913 541-1890
FAX: 913 541-1897
E-mail: orders@heartlandmedspec.com

 Made in the USA
ISO 13485 and FDA Registered

